



## TR 4000 Bariatric Shower Trolley

---

### The TR 4000 Bariatric Shower Trolley Multi Tasking - Endless Possibilities

The TR 4000 Multi Tasking Shower Trolley is designed for lots of power and loads of space for your most difficult patients. The trolley is used for shower, transport, as a changing or treatment trolley, the possibilities are endless.

The patient can be transferred from bed to shower and back again in a safe and secure manner. The height of the base allows mobile patient sling lift transfers. Rechargeable battery powered dual actuators raise and lower the stretcher easily and quietly, with an easy to use hand held remote control. The hand control operates the power tilt in both directions for improved drainage and tilts to a "Trendelenburg" position when desired. An audible tone sounds when the stretcher is level.

The side rails keep the patient in place. Side rails swing down for transfer and lock into up position for patient security for a secure transfer. Locks on each side rail ensure the patients safety. For a tall patient the head end rail can be folded to a horizontal position for extra length. The stretcher has two positions, flat for transfer and transport, tilt for shower for improved drainage. Includes mattress, rechargeable batteries, battery charger, hand held remote control, wedge pillow, drain plug and hose.

---

## TR Bariatric Shower Trolley

### Features

The stretcher is easily and quietly raised by two powerful actuators. The hand held control operates the power tilt and raises and lowers the trolley.

The stretcher can be easily tilted to the "Trendelenburg" position, to alleviate a person in shock. The adjustable stretcher enables a comfortable working height for the

attendant. An audible tone sounds when stretcher is level.

When transporting through narrow doorways the side rails can be pressed together reducing the stretcher width.

The stretcher has rubber bumper wheels on all four corners to protect walls and doorways. Stretcher comes standard with a soft flame retardant

mattress with drainage at the foot end and includes a drain hose.

The trolley is made of durable material, stainless steel and coated steel, making cleaning simple and efficient. Mattress is easily removable for cleaning.

TR reserves the right to make design changes and/or improvements without notice.



Side rails swing down or lock in horizontal position for transfer and up position for transport.



Head end rail folds to horizontal position for extra length.



Convenient single foot control for straight steering and brakes.

## MODEL: TR 4000 Bariatric Shower Trolley

### SPECIFICATIONS:

- Maximum Load 1000 lbs
- Operating range 21" - 33" floor to stretcher
- Overload protection on circuit board
- Battery rechargeable, 24 V
- Construction stainless and coated steel
- Castors central foot control for straight steering and locking brakes
- Hand held remote control raising, lowering, tilt
- Twin actuators battery powered
- Mattress, pillow flame retardant material

DISTRIBUTOR



**IN CANADA:** 75 Dymont Road  
Barrie, Ontario  
L4N 3H6

**Phone:** 1-800-565-7075  
**Fax:** 705-733-3432

**IN U.S.:** 11 N Summerlin Ave, Ste 285  
Orlando, FL 32801  
407-839-0744

**E-MAIL:** [thesales@themedical.com](mailto:thesales@themedical.com)  
**WEBSITE:** [www.themedical.com](http://www.themedical.com)